



# Touchless Urine Sampling Kit

Do you get the right test result?

Up to 42% of urine specimens are contaminated which leads to incorrect results. Sterile and proper urine collection can reduce this by up to 77% (LaRocco et al. 2016).

IVUS Hygienic Urine Sampling Kit are developed to simplify and facilitate clean, hygienic and sterilized urine collection for regular and advanced testing.

Touchless Urine Sampling means, no spill and no mess. Vacuum sampling from the beaker to the tube, secure that the clinical employees don't have to handle the urine sample. The urine will stay in a closed system until testing. The patient receives a combined kit, including a moist skin cleaning towel and a cleaning instruction of urethral. The cleaning with the moist skin towel avoids sample contamination and gives less false positive results. Produced according to CE-IVD guidelines.

- Sterilized, reliable and ready to use for microbiological urine tests
- ✓ Safe, hygienic and easy-to-use
- ✓ Touchless vacuum sampling



#### **User friendly - No contamination**

IVUS kits comes in a box or in a peel-pack-bag and contains all that the patient needs for a single urine sample collection. The kits are easy to open for the patiens. The manual drawings are simple to understand and thereby reduce the risk of contamination.

### Working advantage - Easy & clean

The IVUS urine sampling comes in one kit and is easy to hand out for the clinical worker. The unique design ensures that the urine will stay in a closed system until testing, which minimize spill and mess.





## Hygienic Urine Sampling - Selection

The IVUS kits are packed in peel-pack bags or boxes and contain all you'll need for single-use sample colletion for biological analysis. IVUS is dedicated to deliver the highest standard within urine kits to hospitals and laboratories to provide the best test results, and therefore all kits are sterilized. They are lot traceable and produced according to CE-IVD guidelines.

### **IVUS Urine Kits**

IV-1020	Box kit: urine beaker, 9 ml vacuum U-tube, wet wipe, manual	100x1 set
IV-1021	Bag kit: urine beaker, 9 ml vacuum U-tube with boric acid, wet wipe, manual	100x1 set
IV-1022	Bag kit: urine beaker, 9 ml vacuum V-tube, wet wipe, manual	100x1 set
IV-1023	Bag kit: urine beaker, 9 ml vacuum V-tube with boric acid, wet wipe, manual	100x1 set
IV-1001	Bag kit: urine cup, 10 ml labeled tube, wet wipe, manual	50x1 set
IV-1003	Bag kit: urine cup, 10 ml tube, wet wipe, manual	50x1 set
IV-1005	Bag kit: urine cup - individually paked	75x1 set

### **IVUS Sample Tubes**

IV-2001	10 ml V-tube with cap & label, individually packed - peel-pack bag	200x1 pcs
IV-2002	10 ml V-tube with cap, individually packed - peel-pack bag	100x1 pcs
IV-2003	10 ml V-tube with cap & label, individually packed in two bags - peel-pack bag	50x1 set
IV-2004	10 ml V-tube with cap - 5 pcs packed in two bags - peel-pack bag	35x5 set
IV-2005	10 ml V-tube with cap - 75 pcs packed i one bag	6X75 pcs

### **IVUS Urine Containers**

IV-1040	24h Urine collection container 2000 ml, graduated	52 pcs
IV-1041	24h Urine collection container 3000 ml, graduated, for vacuum tubes	30 pcs
IV-1042	24h Urine collection container 3000 ml. graduated, for vacuum tubes, with 185mm probe	30 pcs

