

EXpell PCR

PCR TUBES & PLATES

CAPP EXpell PCR line comprises of 0.1mL and 0.2mL tubes with plug caps, strips of 8 tubes, strips, plates and optical adhesives. The V-shape bottom allows total sample recovery and the thin-walled wells of polypropylene allow a very good transmission of temperature and minimize condensation.




5100200C



5100210C



5100211C



5100212C



5100213C



5100300C



5100310C



5100311C



5100313C



5100400C

Ordering Information

For standard PCR / qPCR thermocyclers that require regular profile 0.2mL tubes

Cat. No. Description

| | |
|----------|--|
| 5100200C | ExpellPlus PCR tubes 0.2mL, low-binding, bag, 20x1000 pcs. |
| 5100201C | Expell PCR tubes 0.2mL, bag, 20x1000 pcs. |
| 5100210C | Expell PCR strips of 8 tubes 0.2mL, bag, 20x120 pcs. |
| 5100211C | Expell PCR strips of 8 caps (flat), for 0.1mL and 0.2mL tubes, bag, 30x120pcs. |
| 5100212C | Expell PCR strips of 8 tubes w/ single attached caps 0.2mL, 8x120 pcs. |
| 5100213C | Expell PCR plate of 96 x 0.2mL, 8x 25 pcs., non-skirted |
| 5100214C | Expell PCR plate of 96 x 0.2mL, 8x 25 pcs., frosted, sub-skirted (for ABI cyclers) |
| 5100216C | Expell PCR plate of 96 x 0.1/0.2mL, 8x 25 pcs., clear, fully skirted |
| 5100450C | Expell PCR Shell frame grid 96 x 0.2mL for ABI cycler, 8x24 pcs. |
| 5100451C | Rigid shell frame grid 96 x 0.2mL for ABI cycler, 1 pc |

For fast PCR / qPCR thermocyclers that require low profile 0.1mL tubes

Cat. No. Description

| | |
|----------|---|
| 5100300C | Expell PCR tubes 0.1mL, bag, 20x1000 pcs. |
| 5100310C | Expell PCR strips of 8 tubes 0.1mL, bag, 15x120 pcs. |
| 5100311C | Expell PCR strips of 8 tubes w/ single attached caps 0.1mL, 8x120 pcs. |
| 5100313C | Expell PCR plate of 96 x 0.1mL, 8x 25 pcs., non-skirted |
| 5100314C | Expell PCR plate of 96 x 0.1mL, 8x 25 pcs. white semi-skirted (for Roche LightCycler 480 and COBAS Z 480) |
| 5100315C | Expell PCR plate of 96 x 0.1mL, 8x 25 pcs., frosted, sub-skirted (for ABI fast cyclers) |
| 5100400C | Expell PCR Optical Disposable Adhesive, bag of 100 pcs. |